

MOVE WITH CONFIDENCE. LIVE WITHOUT LIMITS.

Innovative Mobility Scooters & Power Wheelchairs
Designed for Life In Motion

Cricket LT Carbon
Model GP309

Weighs ONLY 25.8 lbs.

Without cushion, battery or drink holder.

All Golden scooters and power
wheelchairs are Class II Medical Devices.

 **GOLDEN**[®]
2026 COLLECTION



Enjoy Tight Turns & Total Control

The Golden Buzzaround CarryOn® Power Wheelchair is the **LIGHTEST** Center Wheel Drive Power Chair with the **TIGHTEST Turning Radius** in its class.

It's **ULTRA-COMPACT**,
FOLDS to a **SLIM 10"** with the
BASE WEIGHING ONLY 46 LBS.!
(without seat & battery)



BUZZAROUND CARRYON® FOLDING POWER WHEELCHAIR

Travel-ready power mobility—folds fast, stores small, turns tight.

Folds to 10" | 46 lb. base | Airline-friendly battery

CARRYON
Power
Wheelchair



Scan for video.

TRAVEL-READY



CARRYON® KEY FEATURES:

- Up to 19-mile combined range with 12.5Ah + **FREE 6.5Ah EverRun™ lithium-ion batteries** (travel case included)
- Tight 21" turning radius for easy indoor maneuvering
- Folds to 10" for compact storage (without seat)
- 46 lb. base for easier lifting/transport (without seat & battery)
- 300 lb. weight capacity
- Center-wheel drive for confident control
- Independent 4-wheel suspension and springless articulating chassis for a smoother ride
- **PowerLink™** seat connection for quick setup and breakdown
- Front and rear lights plus a USB port
- Wider, width-adjustable armrests (2") and removable storage basket
- Angle-adjustable footplate with travel handle



Tight 21" Turning Radius

(Easier navigation in tight indoor spaces.)

GP130 Color Shrouds (Comes with red installed and blue included.)

Velvet Red
(Standard)

Coastal Blue
(Standard)

Crisp White
(Optional*)

Canary Yellow
(Optional*)

Stone Grey
(Optional*)

Desert Tan
(Optional*)

Orchid Pink
(Optional*)

Army Green
(Matte)
(Optional*)

Optional Travel Bag Set*



*Upcharge
applies.

BUZZ
AROUND
CarryOn[®]
HD



400 lb.
Weight Capacity

*Heavy-duty travel scooter that
goes the distance and stores easily.*

THE BUZZAROUND CARRYON[®] SCOOTERS

Travel light and ride strong—easy to fold, easy to store, and built for the distance.

300 lb. or 400 lb. capacity | Up to 22-mile range | Airline-friendly batteries

CARRYON
Scooter



Scan for video.

ALL CARRYON[®] SCOOTER KEY FEATURES:

- Up to 22-mile combined range with 12.5Ah + **FREE 6.5Ah EverRun[™] lithium-ion batteries** (travel case included)
- Easy disassembly plus 4-point carry handles for easy transport
- Front & rear LED lights for visibility
- Infinitely adjustable tiller for a comfortable fit
- USB port, drink holder, and cane holder
- Kickstand for compact storage
- Removable front storage basket
- Rear basket for spare battery storage
- Interchangeable color shrouds

GB120 CARRYON[®] KEY FEATURES:

- **300 lb. weight capacity**
- **12" height when folded** (without seat)
- **4" ground clearance**
- Upgraded seat and flip-up armrests

GB220 CARRYON[®] HD KEY FEATURES:

- **400 lb. weight capacity**
- **14" height when folded** (without seat)
- **4.63" ground clearance**
- Brushless motor for more power
- Large 8.25" drive wheels
- Width adjustable, removable flip-up armrests with soft touch texture
- Full-size padded seat



GB120



GB220



GB120 CarryOn[®]



GB220 CarryOn[®] HD



Rear basket
battery storage



Kickstand

TRAVEL-READY



Optional GB120 & GB220 Travel Bag Set*



GB120 & GB220 Color Shrouds (Comes with red installed and blue included.)

Velvet Red (Standard)	Coastal Blue (Standard)	Crisp White (Optional*)	Canary Yellow (Optional*)	Stone Grey (Optional*)	Desert Tan (Optional*)	Orchid Pink (Optional*)	Army Green (Matte) (Optional*)

*Upcharge
applies.

THE BUZZAROUND LX™

Ride in comfort with full-size stability and true portability—built to take you farther in style.

375 lb. capacity | Up to 18-mile range | Easy disassembly

BUZZAROUND LX™ KEY FEATURES:

- 375 lb. weight capacity
- Up to 18-mile range with two U1 batteries
- Disassembles in seconds for easy transport
- Comfort spring suspension for a smoother ride
- Quilted high-back seat with armrests (foldable)
- Infinitely adjustable tiller with EZ Reach adjustment lever
- Ample foot room
- USB charging port & smartphone holder

DISASSEMBLES IN SECONDS

GB149 4-wheel model shown.

3-wheel - 50 lbs.
4-wheel - 57 lbs.

36 lbs.



GB119 & GB149 Color Shrouds

Velvet Red (Standard)	Coastal Blue (Standard)	Crisp White (Optional*)	Canary Yellow (Optional*)
Stone Grey (Optional*)	Desert Tan (Optional*)	Orchid Pink (Optional*)	Army Green (Matte) (Optional*)

*Upcharge applies.



GB119



GB149



THE BUZZAROUND EX™

Full-size comfort with lightweight portability—easy to disassemble, easy to go.

350 lb. capacity | Up to 18-mile range | 4" ground clearance

BUZZ EX



Scan for video.

BUZZAROUND EX™ KEY FEATURES:

- 350 lb. weight capacity
- Up to 18-mile range per charge with two U1 batteries
- 4" ground clearance
- Comfort spring suspension for a smoother ride
- Ultrabright LED headlight & taillight
- Smartphone holder



GB118

BUZZAROUND LX & EX OPTIONAL ACCESSORIES



Y-Cable Charging Harness

Charges both U1 batteries at once.



Battery Docking Station

Off-board charging for U1 batteries.
(Great for travel.)



Bumper Conversion Kit

Provides protection & enables easy vertical storage.

GB118 & GB148 Color Shrouds

Velvet Red (Standard)	Coastal Blue (Standard)	Crisp White (Optional*)	Canary Yellow (Optional*)
Stone Grey (Optional*)	Desert Tan (Optional*)	Orchid Pink (Optional*)	Army Green (Matte) (Optional*)

*Upcharge applies.



GB148

THE BUZZAROUND XL™

Compact travel scooters that disassemble fast and turn tight—built for confident everyday mobility.

300 lb. or 325 lb. capacity | Up to 10-mile range | MicroTurn™ 39" turning radius (XL 4-W)

BUZZ XLSDH



Scan for video.

BUZZAROUND XL™ 3-W KEY FEATURES:

- 300 lb. weight capacity
- Up to 10-mile range with 12Ah batteries
- Easy disassembly for transport
- 3.5" ground clearance
- Angle-adjustable front LED headlight
- Stadium-style contoured seat
- Easy-grip tiller adjustment knob
- USB port on wraparound delta tiller
- Smartphone & drink holders
- Front basket
- Red (installed) & Blue (included) color shrouds + six optional colors

NEW MODEL - GB150

BUZZAROUND XL™ 4-W KEY FEATURES:

- 325 lb. weight capacity
- Up to 8-mile range with 12Ah batteries
- **MicroTurn™ Steering Technology (39" turning radius)**
- 3.5" ground clearance
- Front & rear suspension
- Front & rear LED lighting
- Full-color dash with speed & battery indicators
- Stadium-style mesh seat
- Width-adjustable, flip-up armrests
- USB ports on the tiller
- Smartphone & drink holders
- Front basket
- Red (installed) & Blue (included) color shrouds



microTURN™
STEERING TECHNOLOGY

GB121-STD EXTENDED RANGE

OPTION: 22Ah XR Battery*
(in addition to standard 12Ah)

Packaged separately as an accessory; ships with unit.

*Upcharge applies.

Color Shrouds (No optional shrouds on GB150 & GB152).

Velvet Red (Standard)	Coastal Blue (Standard)	Crisp White (Optional*)	Canary Yellow (Optional*)	Stone Grey (Optional*)	Desert Tan (Optional*)	Orchid Pink (Optional*)	Army Green (Matte) (Optional*)
---------------------------------	-----------------------------------	-----------------------------------	-------------------------------------	----------------------------------	----------------------------------	-----------------------------------	---

THE BUZZAROUND XLSHD™ & XL+™

BUZZ XLSHD



Scan for video.

Heavy-duty support meets compact portability—with premium XL+ upgrades for tighter turns and more control.

325 lb. capacity | Up to 9-mile range (XLSHD) | MicroTurn™ 39" turning radius (XL+)

BUZZAROUND XLSHD™ KEY FEATURES:

- 325 lb. weight capacity
- Up to 9-mile range with 12Ah batteries
- 4" ground clearance
- Two 12Ah batteries (standard)
- Quick & easy disassembly for transport
- Adjustable comfort spring suspension with manually adjustable rear suspension
- Angle-adjustable front LED headlight
- USB port on the wraparound delta tiller
- Smartphone & drink holders
- Front basket

NEW MODEL - GB152

BUZZAROUND XL+™ KEY FEATURES:

- 325 lb. weight capacity
- Up to 9-mile range with 12Ah EverRun™ Lithium-ion battery
- **MicroTurn™ Steering Technology** (39" turning radius)
- 3.5" ground clearance
- Front & rear suspension & LED lighting
- EZ-Connect magnetic charging port and cable
- Stadium-style mesh seat
- Width-adjustable, flip-up armrests
- USB ports on the tiller
- Smartphone & drink holders
- Front basket

SPECIAL GB152 XL+ FEATURES



Color LED Console
with battery meter
& speedometer.



**LED Side Accent
Lighting**
60+ customizable
color options.



EZ Pivot Tiller
Two levers
adjust tiller by
hand or foot.



GB121-SHZ

**BUZZAROUND
XLSHD**

**BUZZAROUND
XL+**

microTURN
STEERING TECHNOLOGY



GB152

39" Turning Radius

EXTENDED RANGE OPTIONS:

GB121-SHZ 22Ah XR Battery*

GB152 24Ah XR Battery*

(in addition to standard batteries)

Packaged separately as an accessory; ships with unit.

*Upcharge applies.

THE BUZZAROUND LT™

Smallest. Lightest. Most portable Buzzaround yet—now with upgraded style and MicroTurn™ maneuverability.

300 lb. capacity | Up to 7-mile range | Disassembles fast | MicroTurn™ (4-wheel)

BUZZAROUND LT™ KEY FEATURES:

- 300 lb. weight capacity
- Up to 7-mile range per charge
- **MicroTurn™ Steering Technology** (GB108)
Turning Radius: 41.75"
- Turning radius (3-wheel): 39"
- Two 12Ah batteries (standard)
- Quick, easy disassembly
- Height-adjustable swivel seat with flip-up, width-adjustable armrests
- Front-mounted, index finger throttle
- Tiller-mounted charging port

BUZZAROUND LT DISASSEMBLED

GB108 4-wheel model shown.



GB107 & GB108 Scratch Resistant Color Shrouds



BUZZAROUND LT™



microTURN™
STEERING TECHNOLOGY

THE COMPANION™

Full-size performance for everyday adventures—stable, maneuverable, and built for comfort.

Up to 500 lb. capacity | Full-size stability | Made in the U.S.A.*

COMPANION



Scan for video.

COMPANION™ KEY FEATURES BY MODEL:

GC340 & GC440:

- Front & rear suspension
- Optional larger battery for increased range

GC340, GC440 & GC540

- Two built-in drink holders on tiller

GC540 (Companion HD):

- Adjustable rear suspension



GC240
350 lb. Capacity



GC340
400 lb. Capacity



GC440
400 lb. Capacity



GC540
500 lb. capacity



Red or Blue or Grey - Customer's Choice



Crimson Red



Azure Blue



Galactic Grey
GC340, GC440
& GC540

**GOLDEN
Companion™
HD**

THE AVENGER™

Rugged, heavy-duty outdoor power with premium comfort.
500 lb. capacity | Up to 18-mile range | 4" ground clearance

AVENGER™ KEY FEATURES:

- 500 lb. weight capacity
- Up to 18-mile operating range per charge
- 4" ground clearance
- 13" oversized, solid front and rear tires
- Adjustable tiller for riding comfort
- Deluxe, reclining high-back captain's seat
- Adjustable front headlight & four-way flashers
- Large front cargo basket (standard)
- Two NF22 batteries (sold separately)



GA541 Colors - Customers Choice



Red



Silver



golden
Avenger™



GA541

THE EAGLE™ & EAGLE ALL TERRAIN

Two outdoor-ready scooters built for confident performance—Eagle for smooth, powerful handling; Eagle All Terrain for off-road capability, premium lighting, and extended-range freedom.

500 lb. capacity | Powerful outdoor performance | All-terrain & extended range options available

BOTH EAGLES KEY FEATURES:

Standard on both Eagle models.

- 500 lb. weight capacity
- Full LED lighting package
- Lockable, water-resistant rear storage compartment
- USB charging port
- Smartphone holder
- Bluetooth® speakers
- Rearview mirrors
- Large front basket
- High-back Captain's seat (comfort-first design)
- Gas-pneumatic adjustable tiller

EAGLE™ KEY FEATURES:

Outdoor stability + comfort tuning for everyday terrain.

- Up to 23-mile operating range per charge
- Powerful 140 Amp controller for confident outdoor performance
- 5" ground clearance
- Two NF-22 batteries (standard)
- Optional 75Ah batteries (up to 32-mile range)

NEW MODEL - GR596

EAGLE AT™ KEY FEATURES:

Suspension + premium controls for off-road confidence.

- Up to 34-mile operating range with two 75Ah batteries
- High-performance 200 Amp controller for exceptional climbing ability
- Front & rear adjustable suspension
- 13" pneumatic turf tread tires
- Full LED lighting upgrade with brake light & turn signals
- Reclining, high-back Captain's seat
- Width, angle & height adjustable armrests
- LCD color console screen with three customizable gauge layouts
- Turning speed reduction sensor
- NFC key tag & conventional key ignition
- Convenient drink holders
- Optional 100Ah batteries

GOLDEN
EAGLE



GR595

GOLDEN
EAGLE
ALL TERRAIN



GR596



GR595 Colors - Customer's Choice
GR596 Colors - Red & Grey Only



USB Charging Port
& Volume Control



Storage Compartment



Eagle Dashboard



EagleAT Dashboard & NFC Key Tag

THE LITERIDER®

Rugged tubular design, easy to clean—and easy to take anywhere.

300 lb. capacity | Lightweight disassembly | Durable tubular frame

LITERIDER® KEY FEATURES:

- 300 lb. weight capacity
- Completely wireless disassembly for easy transport
- Durable tubular metal frame for rugged everyday use
- Up to 12-mile operating range on two 22Ah batteries
- Wraparound delta tiller
- Easy-to-read control panel with backlit battery gauge
- High-back stadium-style seat
- Adjustable armrests
- Two storage baskets
- 9" front & rear tires
- Ultrabright LED adjustable headlight
- Optional 20"×18" stadium-style seat

LITERIDER DISASSEMBLED:

GL141 4-wheel model shown.

Heaviest Piece (Front Half):
3-wheel - 44 lbs.
4-wheel - 52 lbs.



Red or Blue - Customer's Choice



LiteRider[®] GOLDEN



GL111



GL141

THE LITERIDER ENVY®

The easiest-to-transport power wheelchair—lightweight, compact, and designed for life on the go.

300 lb. capacity | Ultra-portable design | Tight indoor turning

ENVY



Scan for video.

LITERIDER ENVY® KEY FEATURES:

- 300 lb. weight capacity
- Lightweight design with easy disassembly for transport
- 28.5" turning radius for excellent indoor maneuverability
- Up to 15.5-mile operating range on 22Ah battery
- Dynamic LiNX® controller
- Large adjustable footplate
- NEW high-back stadium seat with quilted diamond stitch
- Convenient off-board charger
- Large underseat storage basket

LITERIDER ENVY DISASSEMBLED:



GP162 Colors - Customer's Choice



THE STRIDE™

Lightweight, foldable power wheelchair designed for travel, tight spaces, and effortless transport.

300 lb. capacity | Lightweight folding design | Airline-friendly battery

STRIDE



Scan for video.

TRAVEL-READY



STRIDE KEY FEATURES:

- 300 lb. weight capacity
- Lightweight design – weighs as little as 52.6 lbs.*
- Compact folded size for easy storage and transport
- Durable aluminum frame
- Up to 9.3-mile operating range with 10Ah lithium-ion battery
- 23.6" turning radius for tight spaces
- Comfort spring suspension
- Convenient under-seat storage & back seat pocket
- USB charging port
- Stands on end for compact storage

**Excludes removable battery, joystick, cushion, and bag.*

COMPACT & FOLDABLE

Folded dimensions: 42.5" x 24.4" x 15.7"

(Stands upright when folded.)




GP301



*Folded with
Kickstand*

GP301 Color



Black

THE CRICKET™ & CRICKET™ LT

Two ultra-lightweight carbon-fiber power wheelchairs designed for effortless travel and everyday independence.

300 lb. capacity | 100% carbon-fiber frame | Airline-friendly battery | Memory foam cushion

ALL CRICKET™ POWER WHEELCHAIR KEY FEATURES:

- 100% carbon fiber frame
- 300 lb. weight capacity
- Joystick controller with USB charging port
- Contoured memory foam cushion
- Foldable base cushion
- Flip-up armrests
- 35-amp controller
- Off-board charging
- EverRun™ airline-approved lithium-ion battery
- Storage pocket
- Stands on end for easy transport and storage

GP302 CRICKET™ KEY FEATURES:

- Weighs as little as 33 lbs.*
- Up to 9.3-mile operating range with 12Ah battery
- 22" turning radius
- Folded frame lock
- Protective travel bag

*Excludes battery, joystick, cushion, and bag

GP309 CRICKET™ LT KEY FEATURES:

- Weighs as 25.8 lbs.**
- Up to 9.5-mile operating range with 10Ah battery
- 29.5" turning radius
- Spare Battery Mount (Optional battery sold separately)
- Battery Carry Bag & Drink Holder

**Excludes battery, cushion, and drink holder.

GP302 & GP309 Color



Golden
Cricket
CARBON



GP302
Travel Bag



GP302 Folded
with Kickstand



Golden
CricketLT
CARBON



GP309 Battery
Carry Case



GP309 Side Fold
Front View
with Base Cushion



CRICKET



Scan for video.

TRAVEL-READY



CRICKET LT



Scan for video.

MEET THE GOLDEN ALLY™


Compact, foldable power wheelchairs with app-enabled remote control and travel-ready portability.

Smart control app | Tight turning | Foldable portability

OPTIONAL ACCESSORIES*:



Headrest Extension



Cup/Phone Holder



Heated Seat Control



Wireless Remote



Three Ally models—each built for easy transport, confident maneuverability, and smart, everyday independence!

ALL ALLY™ MODEL KEY FEATURES:

- App-enabled remote control (mobile app; optional wireless remote*)
- Airline-friendly EverRun™ lithium-ion battery
- Tight turning for indoor maneuverability
- System monitoring voice alerts
- Spacious leg & foot room
- Ergonomic seating
- LED lighting package
- Folds easily and stands on end for compact storage

TRAVEL-READY



Visit App Store
for Remote Control



*Upcharge applies.

CHOOSE YOUR ALLY™

GP303: lightest | GP304 Pro: most capable | GP305 Sport: performance styling

ALLY



Scan for video.




GP303

GP303 ALLY™ KEY FEATURES:

- 330 lb. capacity
- 28.5" turning radius
- Up to 9-mile operating range
- 49.5 lb. weight (without battery)
- Easy-to-reach underseat storage compartment
- Removable memory foam seat cushion
- Faux leather stitched flip-up armrests
- Water resistant backpack
- USB A & C charging ports

GP303 Ally Colors - Customer's Choice



GP304

GP304 ALLY PRO™ KEY FEATURES:

- 330 lb. capacity
- 37" turning radius
- Up to 13-mile operating range
- 69.6 lb. weight (without battery)
- Easy-to-reach underseat storage compartment
- Omni-directional wheels with self-cleaning V-groove
- Front suspension
- Adjustable width armrests
- Adjustable recline back
- Adjustable leg length footplate
- Removable memory foam seat cushion
- Padded flip-up armrests
- Water-resistant backpack
- USB A & C charging ports

GP304 Ally Pro Colors - Customer's Choice




GP305

GP305 ALLY SPORT™ KEY FEATURES:

- 265 lb. capacity
- 30" turning radius
- Up to 9.5-mile operating range
- 87.5 lb. weight (without battery)
- Omni-directional wheels with self-cleaning V-groove
- Memory foam seating
- Side & rear lighting
- Front suspension
- Power folding
- High-tech styling

GP305 Ally Sport Colors - Customer's Choice



THE COMPASS™

EXPERIENCE CONSTANT CONTACT

Six wheels stay grounded for unmatched stability—designed to reduce rocking and improve curb and threshold performance.

Constant wheel-to-ground contact | Superior stability | Curb & threshold confidence

COMPASS SPORT™ KEY FEATURES:

- 300 lb. weight capacity
- Up to 24-mile operating range on U1 batteries
- 24" turning radius
- Springless, articulating chassis
- Dynamic LiNX® controller
- Large footplate
- Non-marking, flat-free tires
- Foldable, adjustable high-back seat
- Adjustable armrests
- EZ Lock compatible
- Scratch-resistant plastic color shrouds
- Optional 50Ah batteries for up to 34-mile operating range*

golden
**COMPASS™
SPORT**



GP605

COMPASS HD™ KEY FEATURES:

- 450 lb. weight capacity
- Up to 12-mile operating range on NF22 batteries
- 20" turning radius
- Springless, articulating chassis
- Dynamic LiNX® controller
- 75 Amp controller
- Four-post seat mount with slide seat bracket (adjustable seat pitch)
- Non-marking, flat-free tires
- Foldable, adjustable high-back seat
- Adjustable armrests

golden
**COMPASS™
HD**



GP620

Red or Blue - Customers Choice



Red



Blue

*Upcharge applies.

THE BUZZABOUT™

Premium comfort and compact maneuverability—built for everyday travel with distinctive safety lighting.

300 lb. capacity | Tight indoor turning | Motion-sensing LED lighting

BUZZABOUT™ KEY FEATURES:

- 300 lb. weight capacity
- 28.5" turning radius for easy indoor maneuverability
- Up to 11-mile operating range on 22Ah batteries
- High-back stadium seat with contoured support
- Motion-sensing LED running lights on footplate and rear handle
- Quick, easy disassembly into four pieces for transport
- Off-board charging capability
- Dynamic LiNX® controller



**Motion-Sensing LED Lighting
(Front + Rear)**

GP164 Color



GP164

A MESSAGE FROM RICH GOLDEN, CEO



*Rich Golden pictured with his children, third generation Golden Team Members!
Pictured left to right: Rich Jr., Alyssa, Rich & Mike.*

Golden has always stood for comfort, quality, and innovation. Since 1985, Golden has built mobility products that help people move with confidence—designed for dependable performance, lasting comfort, and everyday freedom.

Our scooters and power wheelchairs are engineered for smooth maneuverability and reliable handling, so you can go where you want to go—indoors and out—with greater independence.

Thank you for trusting Golden. We're proud to build America's favorite scooters and power wheelchairs—for your family, and ours.

A handwritten signature in black ink, reading "Rich Golden".

OUR VISION...

To build the most innovative and comfortable products that put you and your family's lives in motion.

TRAVEL LIGHTER

Our best-selling TRAVEL CASE accessory – protective, sleek, and ready to go.

- Stylish, curved hard shell
- Durable, protective construction
- Interior mesh partition
- Adjustable buckle closure
- Size: 13.77" W × 3.14" H × 9.84" D
- Fits any Golden scooter and power wheelchair



MOBILITY ACCESSORIES

Seat Belts Are Available in 2 Sizes for Added Safety!

60" or 78" Seat Belts fit all Golden Scooter & Power Wheelchair Models



Cane Holder




GB120 & GB220
Cane Holder



Safety Flag

(Included with most models.)



Y Cable Charging Harness



Crutch Holder



Scooter & Power
Wheelchair Cover



Oxygen Tank Holder



Walker Holder



Drink Holder



Quad Cane Holder



Rear Basket



Power Elevating Seat*



Travel Case



Smartphone Holder



Pack'n'Go Armbag



Rear View Mirror

*Only available on Companion models.

SPECIFICATIONS

Item	Description	Length	Width	Front Tire	Rear Tire	Center Tire	Casters	Weight Capacity (lbs.)	Maximum Speed (mph)	Ground Clearance	Operating Range ¹ (miles)	Turning Radius	Batteries (2 unless noted)	Weight w/ Batteries	Front Half Weight	Rear Half Weight	Seat Weight w/ Arms	Heaviest Piece	Base Weight
GB119	Buzz LX (3-wheel)	45"	22"	9"	9"	N/A	N/A	375	5	4"	18	40.5"	U1	168	50	36	32	50	N/A
GB149	Buzz LX (4-wheel)	44"	22"	9"	9"	N/A	N/A	375	5	4"	18	54"	U1	174	57	36	32	57	N/A
GB118	Buzz EX (3-wheel)	45"	22"	9"	9"	N/A	N/A	350	5	4"	18	42"	U1	154	46	36	23	46	N/A
GB148	Buzz EX (4-wheel)	44"	22"	9"	9"	N/A	N/A	350	5	4"	18	56"	U1	161	53	36	23	53	N/A
GB121-STD	Buzz XL (3-wheel)	41"	20"	7.5"	8"	N/A	N/A	300	4	3.5"	10	37.5"	12Ah*	93	32	24	21	32	N/A
GB121-SHZ	Buzz XLSHD (3-wheel)	41"	20"	7.5"	8"	N/A	N/A	325	4	4"	9	36.5"	12Ah*	100	35	28	21	35	N/A
GB150	Buzz XL (4-wheel)	42"	22"	7.5"	8"	N/A	N/A	325	4.3	3.5"	8	39"	12Ah	110	42.5	28.2	16	42.5	N/A
GB152	Buzz XL+ (4-wheel)	42"	22"	8"	8"	N/A	N/A	325	4.7	3.5"	9.2	39"	12Ah*	97	42.5	28.2	16	42.5	N/A
GB107	Buzz LT (3-wheel)	41.75"	19-21"	8"	8"	N/A	N/A	300	4	2.5"	7	39"	12Ah	85.5	28	20.5	19	28	N/A
GB108	Buzz LT (4-wheel)	44"	19-21"	8"	8"	N/A	N/A	300	4	2.5"	7	41.75"	12Ah	95	37.5	20.5	19	37.5	N/A
GB120	Buzz CarryOn Scooter	38"	19.5"	6"	7"	N/A	N/A	300	3.6	4"	18	48.75"	12.5+6.5Ah	65.5	N/A	N/A	12	47.5	47.5
GB220	Buzz CarryOn HD Scooter	40"	23.5"	7.25"	8.25"	N/A	N/A	400	4.3	4.6"	14	48.75"	12.5+6.5Ah	77.5	N/A	N/A	20	52.5	52.5
GP130	Buzz CarryOn Power Wheelchair	38"	24"	N/A	N/A	7.75"	5.5"	300	3.85	2"	12.5	21"	12.5+6.5Ah	74	N/A	N/A	22	46	46
GL111	LiteRider (3-wheel)	41.5"	22"	9"	9"	N/A	N/A	300	5	3.75"	12	39"	22Ah	125	44	28	21	44	N/A
GL141	LiteRider (4-wheel)	43.25"	22"	9"	9"	N/A	N/A	300	5	3.75"	12.5	48"	22Ah	133	52	28	21	52	N/A
GC240	Companion (3-wheel)	42.25"	24"	9"	10.75"	N/A	N/A	350	5	3.75"	15	37"	U1	171	45	48	26	48	N/A
GC340	Companion (3-wheel)	46"	24.5"	9"	10.75"	N/A	N/A	400	5	4.5"	19	45"	U1*	207	57	71	29	71	N/A
GC440	Companion (4-wheel)	48"	24.5"	10.75"	10.75"	N/A	N/A	400	5	4.5"	19	56"	U1*	222	67.5	51	29	67.5	N/A
GC540	Companion HD (3-wheel)	48"	26"	10.25"	10.75"	N/A	N/A	500	6	3.75"	16	51"	NF22	225.5	N/A	N/A	50	101.5	N/A
GA541	Avenger	57.5"	24.5"	13"	13"	N/A	N/A	500	7	4"	18 ²	67"	NF22**	265	68	68	55	68	N/A
GR595	Eagle	56.5"	27"	13"	13"	N/A	N/A	500	8	5"	23	69"	NF22*	302.5	N/A	N/A	49.5	179	179
GR596	Eagle AT	63"	31.7"	13"	13"	N/A	N/A	500	11	4"	34	76"	75Ah*	377	N/A	N/A	54	208	208
GP301	Stride	42.5"	24.4"	8"	12"	N/A	N/A	300	3.7	1.6"	9.3	23.6"	10Ah	58	N/A	N/A	N/A	52.6	N/A
GP302	Cricket	36.6"	22.25"	N/A	8.5"	N/A	7"	300	3.7	2.6"	9.3	22"	12Ah	39.3	N/A	N/A	N/A	33.2	N/A
GP309	Cricket LT	36"	24"	N/A	11"	N/A	6.5"	300	3.7	1.25"	9.5	29.5"	10AH	31.3	N/A	N/A	N/A	25.8	N/A
GP303	Ally	41"	24.5"	N/A	10"	N/A	6"	330	5	2.5"	9	28.5"	12Ah	56	N/A	N/A	N/A	49.5	N/A
GP304	Ally Pro	39"	25"	N/A	10"	N/A	10"	330	5	2.5"	13	37"	12Ah	75	N/A	N/A	N/A	69.6	N/A
GP305	Ally Sport	42"	25"	N/A	10"	N/A	10"	265	4	2.5"	9.5	30"	10Ah	92	N/A	N/A	N/A	87.5	N/A
GP162	LiteRider Envy	36"	22.5"	N/A	9"	N/A	6"	300	3.5	2.5"	15.5	28.5"	22Ah	123	37	36	23	37	N/A
GP164	BuzzAbout	36"	23"	N/A	9"	N/A	6"	300	4	2.5"	11	28.5"	22Ah	124	29	33	34	34	N/A
GP605	Compass Sport	39.5"	24"	N/A	N/A	10"	6"	300	4	3.25"	24	24"	U1*	179	N/A	N/A	40	97	97
GP620	Compass HD	40"	25"	N/A	N/A	10"	6"	450	4	3.5"	21 ²	20"	NF22*	226	N/A	N/A	51	101	101

¹ Battery range at 200 lbs., and will vary due to rider weight, drive surface, terrain and battery type.

² Battery range at maximum weight capacity, and will vary due to rider weight, drive surface, terrain and battery type.

* Optional battery upgrade available.

** Batteries must be purchased separately.

Brakes: All Golden Technologies scooters come standard with Dynamic/Electromechanical braking systems.

All Golden Technologies scooters come standard with a sealed transaxle, 24 volt motor.

FDA Class II Medical Devices designed by Golden Technologies are intended to assist persons with chronic physical disabilities.



Scan to view our brochure as a digital flip book!



www.goldentech.com