



Scooter Line Specifications Spring/Summer 2025



	BUZZAROUND											LITERIDER		COMPANION			HEAVY DUTY			
Model	GB119	GB149	GB118	GB148	GB121-STD	GB124-STD	GB121-SHZ	GB124-SHZ	GB107	GB108	GB120	GL111	GL141	GC240	GC340	GC440	GC540	GA541	GR595	
Model Description	Buzz LX (3 wheel)	Buzz LX (4 wheel)	Buzz EX (3 wheel)	Buzz EX (4 wheel)	Buzz XL (3 wheel)	Buzz XL (4 wheel)	Buzz XLSHD (3 wheel)	Buzz XLSHD (4 wheel)	Buzz LT (3 wheel)	Buzz LT (4 wheel)	Buzz CarryOn	LiteRider (3 wheel)	LiteRider (4 wheel)	Companion (3 wheel)	Companion (4 wheel)	Companion (4 wheel)	Companion (3 wheel)	Avenger	Eagle	
Medicare Code	Not Coded	Not Coded	Not Coded	Not Coded	K0800	K0800	K0801	K0801	K0800	K0800	Not Coded	K0800	K0800	K0807	K0807	K0807	K0808	Not Coded	Not Coded	
Weight Capacity	375 lbs...	375 lbs...	350 lbs...	350 lbs...	300 lbs...	300 lbs...	325 lbs...	325 lbs...	300 lbs...	300 lbs...	300 lbs...	300 lbs...	300 lbs...	350 lbs...	400 lbs...	400 lbs...	500 lbs...	500 lbs...	500 lbs...	
Drive Wheels	Rear	Rear	Rear	Rear	Rear	Rear	Rear	Rear	Rear	Rear	Rear	Rear	Rear	Rear	Rear	Rear	Rear	Rear	Rear	
Maximum Speed	5 mph	5 mph	5 mph	5 mph	4 mph	4 mph	4 mph	4 mph	4 mph	4 mph	3.6 mph	5 mph	5 mph	5.1 mph	5.1 mph	5.1 mph	5.7 mph	7 mph	8 mph	
Operating Range¹	18 miles	18 miles	18 miles	18 miles	10 miles/18 miles	10 miles/18 miles	9 miles/18 miles	9 miles/18 miles	7 miles	7 miles	14 miles	12 miles	12.5 miles	19.9 miles	19 miles / 29 miles	19 miles / 28 miles	26.8 miles	29.6 miles²	23 miles	
Ground Clearance Under Deck	4"	4"	4"	4"	3.5"	3.5"	4"	4"	2.5"	2.5"	4"	3.75"	3.75"	3.75"	4.5"	4.25"	6.25"	4"	5"	
Ground Clearance to transaxle/motor	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	2.5"	2.5"	2.5"	N/A	2.5"	3.75"	N/A	3.5	
Turning Radius	40.5"	54"	40.5"	54"	37.5"	45.5"	37.5"	44"	39"	48.5"	48.75"	39"	48"	37"	45"	54"	51"	67"	69"	
Type Batteries	2-U1	2-U1	2-U1	2-U1	2-12 AH	2-12 AH	2-12 AH	2-12 AH	2-12 AH	2-12 AH	12.5 AH	2-22 AH	2-22 AH	2-U1	2-U1	2-U1	2-NF22	2-NF22	2-NF22 / 2.75 AH	
Optional Batteries	N/A	N/A	N/A	N/A	2-22 AH	2-22 AH	2-22 AH	2-22 AH	N/A	N/A	N/A	N/A	N/A	N/A	2-50 AH	2-50 AH	N/A	N/A	2-75 AH	
Standard Colors	Red, Blue	Red, Blue	Red, Blue	Red, Blue	Red, Blue	Red, Blue	Red, Blue	Red, Blue	Red	Red	Red, Blue	Red, Blue	Red, Blue	Red, Blue	Grey, Red, Blue	Grey, Red, Blue	Grey, Red, Blue	Silver, Red	Grey, Red, Tan	
Optional Color Shroud Kits	Grey, Yellow, White	Grey, Yellow, White	Grey, Yellow, White	Grey, Yellow, White	White, Yellow, Grey, Tan, Pink, Green	White, Yellow, Grey, Tan, Pink, Green	White, Yellow, Grey, Tan, Pink, Green	White, Yellow, Grey, Tan, Pink, Green	N/A	N/A	White, Yellow, Grey, Tan, Pink, Green	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	
Freewheel Mode	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	
Electronic Speed Control	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	
Electro-Mechanical Brakes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	
Charger Controller	Off-board 90 Amp	Off-board 90 Amp	Off-board 90 Amp	Off-board 90 Amp	Off-board 50 Amp	Off-board 50 Amp	Off-board 50 Amp	Off-board 50 Amp	Off-board 50 Amp	Off-board 50 Amp	Off-board 50 Amp	Off Board 90 Amp	Off Board 90 Amp	Off Board 90 Amp	Off Board 90 Amp	Off Board 90 Amp	Off Board 120 Amp	Off Board 140 Amp	Off Board 140 Amp	
Dimensions:																				
Length	45"	44"	45"	44"	41"	40"	41"	40"	41.75"	44"	38"	41.5"	43.25"	42.25"	46"	48"	48"	57.5"	56.5"	
Width	22"	22"	22"	22"	20"	20"	20"	20"	19.21"	19.21"	19.5"	22"	22"	24"	24"	24"	30.5"	24.5"	27"	
Height (ground to top of back of seat)	37.5"-39.5"	37.5"-39.5"	37.5"-39.5"	37.5"-39.5"	33.75"-35.25"	33.75"-35.25"	34.25"-35.75"	34.25"-35.75"	31.75"-33.75"	31.75"-33.75"	34.5"	37"- 39"	37"- 39"	37"- 39"	38"- 40"	38"- 40"	42.5"-45.5"	42.5"- 45.5"	49"- 50"	
Height (ground to top of back of headrest)	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	50"- 56"	57"- 62.5"	
Ground to Top of Seat	21.75"-23.75"	21.75"-23.75"	22.75"-24.75"	22.75"-24.75"	21.75"-23.25"	21.75"-23.25"	22.25"-23.75"	22.25"-23.75"	19"-21"	19"-21"	23.5"	23.5"- 25.5"	23.5"- 25.5"	22.25"- 24.25"	23"- 25"	23"- 25"	24.5"- 27.5"	25"- 28"	27"- 28"	
Ground to Top of Deck	6.75"	6.75"	6.75"	6.75"	5"	5"	5.5"	5.5"	4"	4"	6"	6"	6"	6.25"	6.5"	6.5"	8"	6.5"	7.25"	
Deck to Top of Seat	15.75"-17.75"	15.75"-17.75"	16"-18"	16"-18"	16.5"-18"	16.5"-18"	16.5"-18"	16.5"-18"	15"-17"	15"-17"	17.5"	18"- 20"	18"- 20"	16.25"- 18.25"	16.5"- 18.5"	16.5"- 18.5"	16.5"- 19.5"	18.75"- 21.75"	19.75"- 20.75"	
Number of Seat Height Positions	3	3	3	3	3	3	3	3	3	3	N/A	3	3	3	3	3	4	4	2	
Size of Increments	1"	1"	1"	1"	.75"	.75"	.75"	.75"	.75"	.75"	N/A	1"	1"	1"	1"	1"	1"	1"	1"	
Front Axle to Rear Axle	34.7"	33.8"	34.7"	33.8"	31.5"	30.5"	31.5"	30.5"	33.5"	33.5"	28.7"	30"	30"	30"	34.5"	35.5"	36"	38.5"	40"	
Weight of Unit: (assembled)	168 lbs.	174 lbs.	154 lbs.	161 lbs.	93 lbs.	98 lbs.	100 lbs.	105 lbs.	87 lbs.	91.5 lbs.	65.5 lbs.	125 lbs.	133 lbs.	170.8 lbs.	207 lbs.	222 lbs. (U1)	225.5 lbs.	265 lbs.	302.5 lbs.	
Front Half	50 lbs.	57 lbs.	46 lbs.	53 lbs.	32 lbs.	37 lbs.	35 lbs.	40 lbs.	28 lbs.	34 lbs.	N/A	44 lbs.	52 lbs.	45 lbs.	57 lbs.	67.5 lbs.	N/A	68 lbs.	N/A	
Rear Half	36 lbs.	36 lbs.	36 lbs.	36 lbs.	24 lbs.	24 lbs.	28 lbs.	28 lbs.	20.5 lbs.	20.5 lbs.	N/A	28 lbs.	28 lbs.	48 lbs.	71 lbs.	51 lbs.	N/A	68 lbs.	N/A	
Base	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	47.5 lbs.	N/A	N/A	N/A	N/A	N/A	N/A	N/A	179 lbs.	N/A
Seat w/Arms	32 lbs.	32 lbs.	23 lbs.	23 lbs.	21 lbs.	21 lbs.	21 lbs.	21 lbs.	13.5 lbs.	13.5 lbs.	12 lbs.	21 lbs.	21 lbs.	26 lbs.	29 lbs.	29 lbs.	50 lbs.	55 lbs.	49.5 lbs.	
Batteries	2 packs / 24 lbs. ea.	2 packs / 24 lbs. ea.	2 packs / 24 lbs. ea.	2 packs / 24 lbs. ea.	16 lbs./30 lbs.	16 lbs./30 lbs.	16 lbs./30 lbs.	16 lbs./30 lbs.	2 packs / 17.5 lbs. ea.	2 packs / 17.5 lbs. ea.	12.5 AH 5 lbs.	29 lbs./pack	29 lbs./pack	23 lbs. ea.	23 lbs. ea. (U1) 33 lbs. ea. (50 AH)	23 lbs. ea. (U1) 33 lbs. ea. (50 AH)	37 lbs. ea.	37 lbs. ea.	52 lbs. ea. (NF22) 52 lbs. ea. (75 AH)	
Tires (standard):	Solid, low profile	Solid, low profile	Solid, low profile	Solid, low profile	Solid, low profile	Solid, low profile	Solid, low profile	Solid, low profile	Solid, low profile	Solid, low profile	Solid, low profile	Solid, low profile	Solid, low profile	Solid, low profile	Solid, low profile	Solid, low profile	Flat Free Solid	Pneumatic¹	Pneumatic²	
Tires (optional):	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
Front	9"	9"	9"	9"	7.5"	7.5"	7.5"	7.5"	8"	8"	6"	9"	9"	9"	9"	10.75"	10.75"	13"	13"	
Rear	9"	9"	9"	9"	8"	8"	8"	8"	8"	8"	7"	9"	9"	10.75"	10.75"	10.75"	10.75"	13"	13"	
Casters	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
Rear Anti-Tippers	2"	2"	2"	2"	2"	2"	2"	2"	2"	2"	2"	2"	2"	3"	2"	2"	3"	3"	3"	
Standard Seat:	Highback	Highback	Highback	Highback	Stadium Style	Stadium Style	Stadium Style	Stadium Style	Stadium Style	Stadium Style	Stadium Style	Highback	Highback	Highback	Highback	Highback	Extra Wide High Back	Captain's Style	Captain's Style	
Slide Seat	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	Standard	Standard	Standard	N/A	N/A	Standard
Back Height - no headrest	18"	18"	18"	18"	14.75"	14.75"	14.75"	14.75"	12.75"	13.5"	13.25"	16"	16"	16"	15.25"	15.25"	18"	21"	23.75"	
Back Height with headrest	23"-26"	23"-26"	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	27.5" - 31.5"	30.75" - 36.25"
Width x Depth	18" x 16"	18" x 16"	18" x 16"	18" x 16"	17" x 14.5"	17" x 14.5"	17" x 14.5"	16" x 14"	16" x 14"	16" x 14"	15" x 13.75"	18" x 16"	18" x 16"	18" x 16"	18" x 16"	18" x 16"	22" x 20"	20" x 17"	20" x 18"	
Color	Black Vinyl	Black Vinyl	Black Vinyl	Black Vinyl	Black Vinyl	Black Vinyl	Black Vinyl	Black Vinyl	Black Vinyl	Black Vinyl	Black Vinyl	Black Vinyl	Black Vinyl	Black Vinyl	Black Vinyl	Black Vinyl	Black	Black Vinyl	Black Vinyl	
Optional Seats:	Stadium Style	Stadium Style	Stadium Style	Stadium Style	N/A	N/A	N/A	N/A	N/A	N/A	N/A	Stadium Style	Stadium Style	Stadium Style	Stadium Style	Stadium Style	Captain's Style	N/A	N/A	
Back Height - no headrest	16"	16"	16"	16"	N/A	N/A	N/A	N/A	N/A	N/A	N/A	16"	16"	18"	18"	18"	18"	N/A	N/A	N/A
Back Height with headrest	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
Width x Depth	20" x 18"	20" x 18"	20" x 18"	20" x 18"	N/A	N/A	N/A	N/A	N/A	N/A	N/A	20" x 18"	20" x 18"	18" x 16"	18" x 16"	18" x 16"	20" x 20"	N/A	N/A	
Color	Black Vinyl	Black Vinyl	Black Vinyl	Black Vinyl	N/A	N/A	N/A	N/A	N/A	N/A	N/A	Black Vinyl	Black Vinyl	Black Vinyl	Black Vinyl	Black Vinyl	Black	N/A	N/A	N/A
Optional Seats 2:	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	Stadium Seat	Stadium Seat	Stadium Seat	N/A	N/A	N/A	
Back Height - no headrest	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	15"	15"	15"	N/A	N/A	N/A	
Back Height with headrest	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
Width x Depth	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	20" x 18"	20" x 18"	20" x 18"	N/A	N/A	N/A	
Color	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	Black Vinyl	Black Vinyl	Black Vinyl	N/A	N/A	N/A	
Seat Options:																				
Pneumatic Rise Seat Option	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
Power Elevating Seat Option	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	Yes	Yes	Yes	Yes	N/A	N/A	N/A
Footrest:	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
Height adjustable	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
Angle adjustable	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
Warranty:																				
Frame	Lifetime	Lifetime	Lifetime	Lifetime	Lifetime	Lifetime	Lifetime	Lifetime	Lifetime	Lifetime	Lifetime	Lifetime	Lifetime	Lifetime	Lifetime	Lifetime	Lifetime	Lifetime	Lifetime	Lifetime
Drive Train	2 Years	2 Years	2 Years	2 Years	2 Years	2 Years	2 Years	2 Years	2 Years	2 Years	2 Years	2 Years	2 Years	2 Years	2 Years	2 Years	2 Years	2 Years	2 Years	2 Years
Electronics²	2 Years	2 Years	2 Years	2 Years	2 Years	2 Years	2 Years	2 Years	2 Years	2 Years	2 Years	2 Years	2 Years	2 Years	2 Years	2 Years	2 Years	2 Years	2 Years	2 Years
Batteries	1 Year	1 Year	1 Year	1 Year	1 Year	1 Year	1 Year	1 Year	1 Year	1 Year	1 Year	1 Year	1 Year	1 Year	1 Year	1 Year	1 Year	1 Year	1 Year	1 Year

Power Chair Line Specs Spring/Summer 2025



	FOLDABLE						LITE RIDER	COMPASS		BUZZ ABOUT
Model	GP130	GP301	GP302	GP303	GP304	GP305	GP162	GP605	GP620	GP164
Model Description	Buzz CarryOn	Stride	Cricket	Ally	Ally Pro	Ally Sport	LiteRider Envy	Compass Sport	Compass HD	BuzzAbout
Medicare Code	Not Coded	Not Coded	Not Coded	Not Coded	Not Coded	Not Coded	Not Coded	K0823	K0825	Not Coded
Weight Capacity	300 lbs.	300 lbs.	300 lbs.	330 lbs.	330 lbs.	265 lbs.	300 lbs.	300 lbs.	450 lbs.	300 lbs.
Drive Wheels	Center	Rear	Rear	Rear	Rear	Rear	Rear	Center	Center	Rear
Maximum Speed	3.9 mph	3.7 mph	3.7 mph	4.4 mph	5 mph	4 mph	3.5 mph	4 mph	4 mph	4 mph
Operating Range¹	12.5 miles	9.3 miles	9.3 miles	9.3 miles	13 miles	9.5 miles	15.5 miles	24 miles / 34 miles	12 miles⁴	11 Miles
Ground Clearance	2"	1.6"	2.6"	2.5"	2.5"	2.5"	2.5"	3.25"	3.5"	2.5"
Ground Clearance to transaxle/motor	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
Turning Radius	21.25"	23.6"	22"	28.5"	37"	30"	28.5"	24"	20"	28.5"
Type Batteries	12.5 AH Lithium-Ion	10 AH 24V Lithium-Ion	12 AH 24V Lithium-Ion	12 AH 24V Lithium-Ion	12 AH Lithium-Ion	10 AH Lithium-Ion	2-22 AH	2-U1	2-NF22	2-22 AH
Optional Batteries	N/A	N/A	N/A	N/A	N/A	N/A	N/A	50 AH	N/A	N/A
Standard Colors	Red, Blue	N/A	Black Carbon Fiber	Teal, White, Black	Teal, Grey	White, Carbon Fiber Pattern	Red, Blue, Silver, Green, Yellow, Orange	Red, Blue	Red, Blue	Red, Blue
Optional Color Shroud Kits	Grey, Yellow, White	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
Freewheel Mode	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Electronic Speed Control	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Electro-Mechanical Brakes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Charger	Off Board	Off Board	Off Board	Off Board	Off Board	Off Board	Off Board	Off Board	Off Board	Off Board
Controller	40 Amp	50 Amp	35 Amp	40 Amp	40 Amp	40 Amp	40 Amp	40 Amp	75 Amp	40 Amp Remote
Dimensions:										
Length	38"	42.5"	36.6"	41"	39"	42"	36"	39.5"	40"	36"
Width	24"	24.4"	22.4"	24.5"	25"	25"	22.5"	24"	25"	23" at front casters
Height (ground to top of back of seat)	34.25"	38.6"	36.5"	36"	37"	37"	34" - 35"	40" - 42"	38.25" - 42.75"	37.5" - 39.5"
Height (ground to top of back of headrest)	N/A	N/A	N/A	N/A	N/A	N/A	N/A	47" - 52.5"	47.25" - 53.75"	36.75" - 38.75"
Ground to Top of Seat	21"	22"	19.3"	21.5"	20.5"	22.5"	21.5" - 22.5"	19" - 21"	21.25" - 24.25"	21.5" - 23.5"
Ground to Top of Deck	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
Deck to Top of Seat	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
Number of Seat Height Positions	N/A	N/A	N/A	2	N/A	N/A	2	3	5	3
Size of Increments	N/A	N/A	N/A	1"	1"	1"	1"	1"	0.5"	1"
Front Axle to Rear Axle	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
Weight of Unit: (assembled)	74 lbs.	58 lbs.	39.3 lbs.	56 lbs.	75 lbs.	92 lbs.	123 lbs.	179 lbs.	226 lbs.	124 lbs.
Front Half	N/A	N/A	N/A	N/A	N/A	N/A	37 lbs.	N/A	N/A	29 lbs.
Rear Half	N/A	N/A	N/A	N/A	N/A	N/A	36 lbs.	N/A	N/A	33 lbs.
Base	46 lbs.	N/A	N/A	N/A	N/A	N/A	N/A	97 lbs.	101 lbs.	N/A
Seat w/Arms	22 lbs.	N/A	N/A	N/A	N/A	N/A	23 lbs.	40 lbs.	51 lbs.	34 lbs.
Batteries	6 lbs.	3.8	3.8	5 lbs.	5 lbs.	5 lbs.	27 lbs. pack	22 lbs. ea. / 30 lbs. ea.	37 lbs. ea.	28 lbs.
Tires (standard):	Flat Free	Flat Free	Flat Free	Pneumatic	Pneumatic + Omni Directional	Pneumatic + Omni Directional	Flat Free	Flat Free	Flat Free	Flat Free
Tires (optional):	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
Front	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
Rear	7.75"	12"	8.5"	10" Pneumatic	10" Pneumatic	10" Pneumatic	9"	10"	10"	9"
Casters	5.5"	8"	7"	6" Solid	10" Omni-Directional	10" Omni-Directional	6"	6"	6"	6"
Rear Anti-Tippers	N/A	2"	1.6"	2"	2"	2"	(Rear) 3"	N/A	N/A	2"
Standard Seat:	Stadium Style	Sling Back with Removable Seat Cushion	Mesh Back with Removable Seat Cushion	Stadium Style	Stadium Style	Stadium Style	Highback	Captain's Style	Captain's Style	Highback
Slide Seat	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
Back Height - no headrest	19.3"	16"	19.3"	18"	19"	18"	16"	23"	19"	18"
Back Height with headrest	N/A	N/A	N/A	N/A	N/A	N/A	N/A	30.5" - 34.5"	26.5" - 30"	30.5" - 34.5"
Width x Depth	17.3" x 16.5"	17" x 18"	17.3" x 16.5"	17" x 16"	21" x 17"	17" x 17"	18" x 16"	18" x 18" / 20"	20" x 20" / 18"	18" x 16"
Color	Black	Black	Black	Black Fabric	Black Fabric	Black / Grey Fabric	Black Vinyl	Black Vinyl	Black Vinyl	Black Vinyl
Optional Seats:	N/A	N/A	N/A	N/A	N/A	N/A	Stadium Style	Captain's Style	Captain's Style	Captain's Style
Back Height - no headrest	N/A	N/A	N/A	N/A	N/A	N/A	16"	19"	19"	18"
Back Height with headrest	N/A	N/A	N/A	N/A	N/A	N/A	N/A	26.5" - 30"	26.5" - 30"	30.5" - 34.5"
Width x Depth	N/A	N/A	N/A	N/A	N/A	N/A	20" x 18"	16", 18", 20" or 22" x 16", 18", or 20"	18", 20", 22" or 24" x 18" or 20"	17" x 16"
Color	N/A	N/A	N/A	N/A	N/A	N/A	Black Vinyl	Black Vinyl	Black Vinyl	Black Vinyl
Optional Seats 2:	N/A	N/A	N/A	N/A	N/A	N/A	Premium	N/A	N/A	N/A
Back Height - no headrest	N/A	N/A	N/A	N/A	N/A	N/A	17"	N/A	N/A	N/A
Back Height with headrest	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
Width x Depth	N/A	N/A	N/A	N/A	N/A	N/A	17" x 16"	N/A	N/A	N/A
Color	N/A	N/A	N/A	N/A	N/A	N/A	Black Vinyl	N/A	N/A	N/A
Seat Options:										
Pneumatic Rise Seat Option	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
Power Elevating Seat Option	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
Footrest:	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Height adjustable	N/A	N/A	N/A	Yes	Yes	No	Yes	2.5" - 4.5"	3" - 7"	No
Angle adjustable	N/A	N/A	N/A	Yes	Yes	Yes	Yes	0 - 10 degrees	0 - 35 degrees	Yes
Warranty:										
Frame	5 Years	5 Years	5 Years	5 Years	5 Years	5 Years	5 Years	Lifetime	5 Years	5 Years
Drive Train	13 Months	13 Months	13 Months	13 Months	13 Months	13 Months	13 Months	13 Months	13 Months	13 Months
Electronics²	13 Months	13 Months	13 Months	13 Months	13 Months	13 Months	13 Months	13 Months	13 Months	13 Months
Batteries	12 Months	12 Months	12 Months	12 Months	12 Months	12 Months	12 Months	12 Months	12 Months	12 Months